



CHROMALIVE[®]

Holistic Colour Balancing Concept
for Healthcare, Beauty and Wellbeing

Auramed[®]

CHROMALIVE®

The History of Colour Light Therapy Research

Throughout history man has used the sun's energy as medicine. Healers would bring the sick out into the sun in order to improve their health. The ancient Egyptians, Greeks and Chinese all used colour therapy, although they used different methods of application of the colour energy. The basis for the scientific explanation of colour therapy came from Sir Isaac Newton with the discovery of spectrum analysis. Johann Wolfgang von Goethe examined the effect of colours on the human psyche and proved the close connection between colours and feelings.

In 1903, Professor Niels Finsen received the Nobel Prize in Medicine for his research in the area of colour therapy. Finsen proved that colours are independent energies and that specific colour vibrations cause different reactions in the human body and the

psyche. In 1933, Dr. Dinshah Ghadiali wrote about his findings and medical research on the application of colours in a Spectro-Chrome Metry Encyclopedia. Since then more than 6,000 cases of colour healing have been documented in American clinical trials.

Dr. Dinshah Ghadiali created combinations of colour light treatments for a variety of conditions, such as burns, fractures, asthma, neurodermatitis, diabetes, psychic disturbances, etc.

Many authors continue to cite Ghadiali's research data and many medical practitioners refer to it. Ghadiali greatly valued having an aura diagnosis before treatment, which he determined through the use of various instruments and methods that he developed.



Coloured Light in Medicine, Wellness and Cosmetics

Since ancient times it was known to the therapists that the stimulation of reflex zones, acupuncture points and meridians had a direct influence and harmonizing effect on organ functions and the psyche. Since the cells react to colour stimuli and pass them on to the brain and the autonomic nervous system, the vital function can be activated and the organism as well as the psyche can be addressed to in a holistic way.



The Healing Power of Light

The Nobel prize winner Prof. Dr. Finsen used colour light to heal a variety of dermal and infectious diseases. In traditional medicine, infrared and blue light are also used frequently to treat different ailments. In modern clinics blue light is often used for the treatment of icterus, which is a common disorder in babies. Dentists use blue light for pain relief and prevention of infections. Even the flight industry recommends the use of green light to calm anxious and nervous people and to increase endurance and energy.



Penlight Sets



Penlight Set Elite silver
Item No. 3030 (glossy, silver-chrome coloured)



Penlight Set Elite gold
Item No. 3020 (24 carat gold-plated)



Penlight Set Basic 2
Item No. 3040 (chrome coloured)

Penlight Set Elite silver Item No. 3030 (glossy, chrome coloured)

This package includes:

- 1 Chromalive® Colour Penlight (glossy-chrome) with an ultra-light special LED (light emitting diode), an integrated battery housing with 3 AAA batteries (1.5 V each) and an on/off turn-switch.
- One broad head - glossy chrome, with a laminar quartz crystal (26 mm diameter). This high quality crystal has a special finish created through a medical fine grinding. The solid crystal surface feels soft and velvety on the skin and is also convenient for pressure point applications. This crystal head is an ideal tool for colour light treatment and massage of wider skin areas and meridians.
- One punctiform head - glossy chrome, with a pointed tip made of high quality crystal created through a medical fine grinding. The solid crystal surface has a diameter of 5.5 mm and the tip is slightly rounded, allowing you to apply accurate pressure

Penlight Set Basic 2 Item No. 3040 (chrome coloured)

This package includes:

- 1 Chromalive® Colour Penlight (anodised in aluminium colour) with an ultra-light special LED and two AAA batteries. The Colour Penlight (three parts) consists of a filter head with a pointed tip made of high quality crystal created through a medical fine grinding. The solid crystal surface has a diameter of 5.5 mm and the tip is slightly rounded, allowing you to

on a specific skin area, without injuring the skin surface. Applications for this crystal head include dissolving tensions and hardenings in the reflex zones, lines of the hand as well as colour light treatment of acupuncture points and scars.

- Seven colour filters in the colours red, orange, yellow, green, turquoise, blue and magenta. The filter frame consists of aluminium and the colour inlays are made of a stiff plastic film.
- One carrying case (silver-grey) made of high quality plastic with embossed printing and a fitted polyethylene foam insert.
- One instruction manual.

Penlight Set Elite gold Item No. 3020 (24 carat gold-plated)

- It has the same specifications as the penlight Elite - silver-chrome, but the body and the heads are with a 24 carat gold-plated finish. (Price on request!)

apply accurate pressure on a specific skin area, without injuring the skin surface. Applications for this crystal head include dissolving tensions and hardenings in the reflex zones, lines of the hand as well as colour light treatment of acupuncture points and scars.

- 7 interchangeable colour filters in the following colours: red, orange, yellow, green, turquoise, blue, magenta.
- An attractive carrying case with a foam fitted inlay.
- One instruction manual.

The innovative Colour Therapy Concept by Auramed®

Based on the scientific researches of Dr. Ghadiali Auramed developed a professional concept of aura analysis and colour harmonisation for therapy, wellness and beauty. This is distributed internationally and available in several languages.

The Auramed Colour Therapy Concept consists of the following steps:

1. Aura Measuring:

Establishing an aura measurement of 43 organ zones using the Biopulsar-Reflexograph® biofeedback sensor.

2. Analysis and Coaching:

On the basis of the Aura and Biofeedback Software.

3. Chromalive® Colour Light Harmonisation:


A method of treatment on the palm and/or facial reflex zones with the Chromalive® Software and Penlights specially developed by Auramed®.



Auramed® offers various sensor designs and software modules.

Sensors as well as software modules can be combined according to individual requirements and upgraded accordingly.

For working with the Biopulsar-Reflexograph® user seminars and web training as well as literature are offered.



Colour Therapist Medical



Colour Therapist Medical

Item No. A5040

The ideal basic equipment for the work in the field of Therapy, Wellness, Proactive Healthcare and Beauty.

The package **Colour Therapist Medical** contains:

- **One Biopulsar-Reflexograph® Favorit Biofeedback Sensor** with 2 different palm sizes. The measuring electrodes for 43 organ zones on each hand are 24 carat gold plated and the housing consist of white polyethylen. An attractive hard shell case with the engraving Auramed® and matching foam insert provides the sensor protection and is therefore suitable for mobile use.

Also included is a technical manual, a mousepad with the description of the reflex zones, a user manual and the booklet: *The Energy Language of the Organs*.

- **The Software** is delivered on CD and consists of the analysis modules **Biofeedback, Aura** and the therapy software **Chromalive®**.

The Chromalive® software contains information on each organ in text and images and gives detailed instructions on how to harmonize the facial and hand reflex zones with coloured light, based on the measured biofeedback values. Furthermore, recommendations are given on essential oils and positive affirmations to support the therapy.


- **A Chromalive® Penlight Set Elite silver** with 2 different quartz crystal heads and 7 colour filters in a padded hard case with the inscription Auramed®.

Software

Excerpts from the **Colour Therapist Medical** Software Modules:



Biofeedback



Aura



Facial Zones



Hand Zone